**EN** These instructions for use have been translated and are available at the following link: [www.c-tech-implant.com](http://www.c-tech-implant.com). Otherwise, scan the following QRcode. Alternatively, request a review of your desired IFU at the following email address [quality@c-tech-implant.com](mailto:quality@c-tech-implant.com). For further information, read the submission methods in the "Information for Customers" paragraph.

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**C-TECH DENTAL IMPLANT SYSTEMS**  
**INSTRUCTIONS FOR THE USE OF DENTAL IMPLANTS**

**Product identification and intended use**

EL, CP, ND, GL and BL dental implants are medical devices made for use exclusively in the oral cavity, in patients suffering from total or partial mandibular and/or maxillary edentulism (absence of natural teeth). The devices are inserted into the bone to act as a substitute for the natural tooth root that will provide an anchor for the permanent dental crown to replace missing natural teeth both functionally and aesthetically. They are replaceable devices that come into contact with bone and soft tissue in the mouth. They have a long service life if the patient health conditions are ideal, and all the necessary maintenance interventions are carried out.

The external structure of an implant has a conical shape with a thread that allows it to make its way into harder bone, maintain its stability in softer bone and not damage the more sensitive areas of the mouth such as the canal nerve or Schneider's membrane (maxillary sinus). Implant threads are therefore different depending on the indications for use and the biological structures that patients present. The threads and apices can be more aggressive when they have to be inserted into a harder bone. On the other hand, they may be less aggressive if they have to be inserted into areas of softer bone or near more delicate areas (nerves, membranes).

All implants are equipped with an internal connection system that allows the devices intended for prosthetic rehabilitation to be solidly anchored to the implants, thus allowing the restoration of missing elements. The EL, CP and BL implants have the same connection, which means that the same prosthetic components can be used on both. Before any such use, the compatibility of the parts should be assessed at the level of gingival tissues and for the type of prosthetic rehabilitation to be performed.

EL, CP, BL and GL prosthetic components can be used on all diameters and lengths of the implant lines available. More details on the diameters and lengths available can be found under "Supply, Diameters and Lengths".

Original components and products from C-Tech Implant must be used for prosthetic rehabilitations of the EL, CP, ND, BL and GL lines. C-Tech disclaims any liability for the use of non-original instrumentation. It is recommended that only appropriate instruments supplied by C-Tech Implant, which are available individually or in special surgical organizers, be used for the surgical insertion procedure. C-Tech Implant disclaims any liability for the use of non-original equipment. For more details on the indications for use of each line, please see the section on 'Indications and Contraindications'.

EL and CP implants are manufactured from grade 4 titanium and are surface-treated with acidification. ND, BL and GL implants are manufactured from grade 5 titanium and subjected to sandblasting surface treatment. Cover screws are manufactured in Titanium Grade 5. The types of titanium used for manufacture comply with harmonised standards.

**IDENTIFICATION OF THE MANUFACTURER:** The manufacturer of the EL, CP, ND, BL and GL implant lines is:

**C-TECH IMPLANT SRL**, Via Ravenna, 382 40018 – San Pietro in Casale (BO) – Italia, Tel. +39 051.6661817, email: [info@c-tech-implant.com](mailto:info@c-tech-implant.com), [implant.com](http://implant.com), website: [www.c-tech-implant.com](http://www.c-tech-implant.com)

**SUPPLY, DIAMETERS AND LENGTHS**

Dental implants are supplied in plastic vials inside heat-sealed, sterile blister packs (secondary packaging) with each containing 1 implant and 1 cover screw (H. 0). The blister and outer carton are duly labelled and bear the code, expiry date, batch number and UDI code. The labels on the blister pack are removable and must be used by the medical personnel to be affixed to the patient implant card, which is also provided inside the implant case to then be delivered to the patient. They must also be attached to the patient's medical record. The Instructions for Use are also supplied with the implant and the patient implant card. The implant card must be filled in in all its fields according to what is requested and based on the information held by the medical staff. Please check the "Table 1" for a better understanding of the symbols on the medical device labels and the "Table 2" for a better understanding of the symbols on the patient implant card.

Implants are available in the following diameters and lengths:

LINE	DIAMETERS AND LENGTHS
ND	Ø 3.1 mm L. 9/11/13 mm
EL	Ø 3.8 mm L. 8/9/11/13/15 mm Ø 4.3 mm L. 7/9/11/13/15 mm Ø 5.1 mm L. 7/9/11/13/15 mm Ø 6.0 mm L. 7/8/10/12/14 mm Ø 7.0 mm L. 7/8/10/12/14 mm
CP	Ø 3.8 mm L. 8/9/11/13/15 mm Ø 4.3 mm L. 7/9/11/13/15 mm Ø 5.1 mm L. 7/9/11/13/15 mm Ø 6.0 mm L. 7/8/10 mm
BL	Ø 3.8 mm L. 7/8/9/11/13/15 mm Ø 4.3 mm L. 6/7/8/9/11/13/15 mm Ø 5.1 mm L. 6/7/8/9/11/13/15 mm
GL	Ø 3.5 mm L. 8/9/11/13/15 mm Ø 4.3 mm L. 8/9/11/13/15 mm Ø 5.1 mm L. 8/9/11/13/15 mm

**STORAGE AND USE:** Implants should be stored in a cool, dry place at room temperature, away from heat, sunlight and water. Sterile products must be used no later than the expiry date stated on the labels. Do not use devices from open or damaged packages or blisters whose sterile status can no longer be ensured.

**WARNINGS:** • C-Tech Implant disclaims all liability for damages caused to the patient by improper use or use other than as described in this instruction sheet or use by unqualified and/or untrained personnel. • The law allows the surgical use of dental implants only and exclusively by qualified surgeons and dentists. Use by other persons is absolutely forbidden. • The manufacturer disclaims any liability for the use of unmarked medical devices not supplied by C-Tech Implant. • The surgical techniques required to place dental implants are highly specialised and complex procedures. Professionals will have to attend and complete courses to be ready for implantology techniques. • C-Tech Implant medical devices are provided with information materials and the company organises training courses worldwide. Participation in a specialised course is therefore strongly recommended. • The use of improper techniques can lead to patient injury, implant failure, bone loss and post-operative complications. • Dental implants are supplied in unopened, sterile packages and are ready for use. These devices are disposable, non-reusable and cannot be re-sterilized. Any re-use or re-sterilisation of medical devices is strictly prohibited. This is because the functional performance and safety of the device expected by the manufacturer cannot be guaranteed, exposing the patient to the risk of loss of device performance, injury, pain and/or serious risk of contamination and cross-infection. • Do not use a device in packaging that has been

tampered with. • Do not use the device beyond the expiry date • Due to potential effects of anaesthesia do not use on pregnant women. • Use the devices on adults or on people who have reached the age of majority. • The device has not been evaluated in paediatric/adolescent patients and is not recommended for use on children. Routine treatment is not recommended until the end of the juvenile growth of the maxillary/mandibular bone has been properly documented. • Inform the healthcare professional of the presence of dental implants in case it is necessary to undergo magnetic fields, diagnostic investigations, therapeutic treatments or other procedures. In most cases, titanium dental implants do not interfere with magnetic fields and there are no known contraindications of any kind for therapeutic treatments or diagnostic investigations.

**INDICATION FOR USE**

**ND Line:** ND implants are suitable for supporting single, fixed or removable dentures or bridges, to rehabilitate masticatory function. To achieve this result, a one- or two-stage surgical technique can be used in combination with immediate, early or delayed loading protocols, as long as sufficient primary stability and an appropriate occlusal loading for the selected technique have been acknowledged. **Use:** in D1, D2, or D3 bone to replace upper lateral incisors, lower lateral and central incisors. They may be used on premolars but only after a careful assessment by the doctor and of the prosthetic protocol. In this case, it is not recommended that stumps that are not straight be used. **Do not use:** On upper and lower posterior sectors and in D4 bone.

**EL & CP Line:** EL and CP implants are indicated to support prostheses ranging from single tooth to full arch fixed-removable dental procedures to rehabilitate the masticatory function. To achieve this result, a one- or two-stage surgical technique can be used in combination with immediate, early or delayed loading protocols, as long as sufficient primary stability and an appropriate occlusal loading for the selected technique have been acknowledged. **Use:** Ø 3.8 on all types of bone for teeth from the second premolar to the upper and lower central incisor. For use on the posterior sector, assess according to bone availability. Use Ø4.3 from the first premolar to the second molar. Use Ø5.1 to Ø7.0 on molars and evaluate according to bone dimensions. **Do not use:** on narrow bone ridges.

**BL Line:** BL implants are indicated to support prostheses ranging from single tooth to full arch fixed-removable dental procedures to rehabilitate the masticatory function. To achieve this result, a one or two-stage surgical technique can be used in combination with immediate, early or delayed loading protocols, as long as sufficient primary stability and an appropriate occlusal loading for the selected technique have been acknowledged. **Use:** Ø 3.8 on all types of bone for teeth from the second premolar to the upper and lower central incisor. For use on the posterior sector, assess according to bone availability. Use Ø4.3 from the first premolar to the second molar. Use Ø5.1 on molars and evaluate according to bone dimensions. **Do not use:** on narrow bone ridges.

**GL Line:** GL implants are indicated to support prostheses ranging from single tooth to full arch fixed-removable dental procedures to rehabilitate the masticatory function. To achieve this result, a one- or two-stage surgical technique can be used in combination with immediate, early or delayed loading protocols, as long as sufficient primary stability and an appropriate occlusal loading for the selected technique have been acknowledged. **Use:** Ø 3.8 on all types of bone for teeth from the second premolar to the upper and lower central incisor. For use on the posterior sector, assess according to bone availability. Use Ø4.3 from the first premolar to the second molar. Use Ø5.1 on molars and evaluate according to bone dimensions. **Do not use:** on narrow bone ridges.

**CONTRAINDICATIONS: Absolute:** Established allergies to pure titanium grade 4, to Ti-6AL-4V titanium alloy, recent myocardial infarction, valvular prostheses, severe renal disorders, severe liver disorders, osteomalacia, treatment-resistant diabetes, recent high-dose radiotherapy, severe chronic alcoholism, drug use, terminal illnesses and tumours. **General:** Uncooperative patient. Intraoral contraindications such as: Anatomical impediments, insufficient bone tissue or poor bone quality. **Related:** Chemotherapy, mild renal disorders, mild liver disorders, endocrine disorders, psychological disorders or psychosis, lack of understanding and motivation, AIDS, seropositive cases, prolonged use of corticosteroids, calcium-phosphorus metabolism disorders, erythropoietic disorders.

**WARNING!** The success of an implant cannot be 100% guaranteed if the indications for use of the product are not followed and the appropriate surgical/management procedures are not put in place. Close collaboration among the surgeon, the prosthodontist and the dental laboratory is essential for successful implant treatment. When using a new device or treatment method for the first time, it is advisable to work alongside a colleague who is experienced in the use of the new device/treatment method to try to avoid possible complications. C-Tech Implant will gladly put you in contact with a product manager, its network of consultants and key opinion leaders.

**PATIENT SELECTION:** The medical devices covered by this leaflet are to be used on patients who have been diagnosed with total and/or partial mandibular and/or maxillary edentulism. The device has not been evaluated in paediatric/adolescent patients and is not recommended for use on children. Routine treatment is not recommended until the end of the juvenile growth of the maxillary/mandibular bone has been properly documented. Ethnicity has no impact on the use of these medical devices. Not all edentulous patients are amenable to implant treatment, which is why it is necessary for medical personnel to subject the patient to an appropriate and careful initial selection process and then a precise diagnosis. Patient assessment prior to an implantology surgical procedure is extremely important. This includes assessment of the patient's general health, oral hygiene status, dental and periodontal condition and their anatomical acceptability. Taking the patient's medical and clinical history is essential. Allergies to titanium or to titanium alloys are rare but possible. It is advisable to check in advance with patients not to have allergies of this type. The manufacturer provides the technical data sheets of the materials upon request for a deeper investigation. Particular attention should be paid to patients with local or systemic factors that could interfere with or slow down the healing processes of bone or soft tissue and osseointegration (e.g. smoking, poor oral hygiene, uncontrolled diabetes, orofacial radiotherapy, steroid therapy, or infections in the surrounding bone). Special attention should be paid to patients undergoing bisphosphonate therapy. Patients with a hyperfunction (bruxism) may be at greater risk of implant failure or fracture. In general, for patient selection and prosthetic design the dentist should consider all the biological and physical conditions of the patient. The elements for this analysis are acquired during the doctor's studies and career path.

**PRE-OPERATIVE PLANNING AND TREATMENT:** Preparation for the surgical procedure includes: - General medical, psychological and dental history, general medical examination, clinical (complete blood count) and radiological examinations, CT scan and, if necessary, consultation with the patient's primary care doctor. - Patient information (indications, contraindications, clinical picture, possible undesirable effects, need for periodic follow-ups) - Dental hygiene plan, with any required periodontal procedures. - Evaluation of the procedure and surgical technique based on the treating dentist's studies, knowledge and experience. Since surgical techniques are highly specialised and complex procedures, it is strongly recommended that practitioners complete a training course to learn established placement techniques in oral implantology. Failure to seek proper training or to follow the instructions for use (IFU) can affect the clinical success of the implant, such as implant failure, bone loss or postoperative complications, which could lead to pain for the patient and put his or her health at serious risk. Surgical placement and subsequent reconstruction with dental implants require complex procedures that may involve risks. It is the doctor's responsibility to inform the patient of any risks and treatment alternatives (including a 'no treatment' option) before performing any procedure.

**INDICATIONS FOR USE AND SURGICAL PROCEDURE:** Surgical techniques for implants are generally taught in university undergraduate programmes in dentistry. Nevertheless, the recommendations in this section should be considered. Users are encouraged to attend professional training courses and/or to read the information materials provided by the company and available on the website at [www.c-tech-implant.com](http://www.c-tech-implant.com). At the bottom of this IFU there is a QR code that links to the page on the information material site.

**WARNING!** The surgical procedure for inserting implants into the oral cavity must take place in a field that is as sterile as possible. Medical personnel are to be equipped with all necessary personal protective equipment (e.g., gloves, mask, gown, etc.) during surgery and whenever

**handling any medical devices.** To proceed with the surgical procedure, medical personnel should be equipped with the surgical instrumentation sold by C-Tech Implant, which must be compatible with the implant line and the surgical technique to be used. Any use or reuse of surgical instruments without having carefully and thoroughly read the specific package leaflet cannot ensure the functional and safety performance of the devices provided by the manufacturer. The patient could therefore be exposed to risks of non-mechanical product performance and/or serious contamination. The surgical instruments supplied by the manufacturer must be cleaned, disinfected and sterilised before first use and after each re-use in accordance with the IFU sheet accompanying those instruments. All surgical instruments must be checked for cleanliness and sterility before each use to ensure that they are always ready to be used. Never use any devices with anomalies or wear as this could compromise the patient's health and safety. Inspect all surgical instruments before use. The choice of the surgical procedure to be performed is to be assessed and decided upon by the doctor, in collaboration with the dental technician, after having carefully taken the patient's medical history and only once the prosthetic project has been properly planned and organised. There are two surgical procedures that can be performed with C-Tech devices: traditional surgery and guided surgery (the latter are applicable only for the EL, CP and ND lines).

**Traditional surgery:** Carry out the tissue incision procedure using the most appropriate technique for management of the case and which is best suited to the patient's physical and biological characteristics. The moment when the implant site is milled is crucial for a successful outcome. It is important that each step is followed to the letter: irrigate constantly while drilling and then set an appropriate cutting speed for the instrument being used and for the bone encountered in the mouth. The speed selected by the dentist when preparing the implant site depends on the preparation technique chosen. The wide range available allows the dentist to choose the most appropriate technique whilst remaining within values that do not lead to bone tissue damage. For additional details on the surgical procedure, please attend the courses organised by the company and read the materials available on the website.

When preparing the implant site, it is necessary to start by creating a hole with a locator drill (100-750 RPM), and then proceed using the pilot drill (100-750 RPM), with the desired stop length, to drill the hole for the chosen implant length. The EL & CP implant should be inserted 1 mm below the bone crest or at bone level. So, the final length of the hole should also take this measurement into account. Then continue sequentially using the drills measuring  $\varnothing$  3.1-  $\varnothing$  3.5-  $\varnothing$  4.3-  $\varnothing$  5.1-  $\varnothing$  6.0-  $\varnothing$  7.0 (100-750 RPM) until the diameter of the chosen implant is reached. To facilitate implant insertion, in the case of compact bone (D1, D2), a step with the counterbore (200-600 RPM), with the hard bone drills (100-750 RPM) and also with the bone tap (15-30 RPM) is recommended.

**Guided surgery (only for EL, CP  $\varnothing$  3.8 to  $\varnothing$  5.1 and ND implants):** To successfully use guided surgery, the dentist must have good previous experience in traditional surgery and also very good knowledge of anatomy. For proper patient selection for this type of surgical procedure, follow the instructions given in the chapter on 'patient selection' and apply what you learned during your studies. It is important that each step is followed to the letter: irrigate constantly while drilling and then set an appropriate cutting speed for the instrument being used and for the bone encountered in the mouth. The speed selected by the dentist when preparing the implant site depends on the preparation technique chosen. The wide range available allows the dentist to choose the most appropriate technique whilst remaining within values that do not lead to bone tissue damage. For additional details on the surgical procedure, please attend the courses organised by the company and read the materials available on the website.

Remove the soft tissue portion using a tissue punch (100-800RPM), and then in case of irregular bone ridges, level the bone with the bone profiler (150-400RPM). Then proceed by drilling a pilot hole for the drills with a locator drill (400-800 RPM). Then the dentist should proceed to drill sequentially (400-800 RPM), using as many drills as as many diameters as necessary, until reaching the diameter and burr length that correspond to the chosen implant length. In the case of compact bone (D1, D2), to facilitate implant insertion, a step using hard bone drills (400-800 RPM) is recommended.

**POSITIONING THE IMPLANTS:** Irrigate the entire implant site with sterile water or saline solution. Using gloves, remove the plastic tube from the blister pack and remove the top cap. Do not take the implant out of the tube. Using an implant wrench, either a ratchet or contra-angle wrench, first, engage the implant connection to remove it from the tube without touching it with your hand. Screw the implant down into the bone with the contra-angle torque wrench (between 35 and 40 Ncm). Then if necessary use the torque wrench (between 35 and 40 Ncm) to screw down the implant, clockwise, until it is completely inserted in the bone. If, in order to achieve complete insertion, greater force is required, the implant must be carefully unscrewed and removed. Then, use the hard bone burr once again to facilitate its insertion. A locking screw is included above the implant cap. This screw is to be screwed down tightly on the implant before the flap is sutured (if it is open). A radiographic check of the implant placement is recommended prior to flap closure (if open).

**POST-OPERATIVE RECOVERY:** The patient should be supervised and instructed to follow a post-operative course appropriate to the surgery undergone. This includes the application of cold compresses, and no physical exertion and/or sauna for the first 24 hours. Alcohol, nicotine, tea or coffee should be limited for the first few days, and the patient should avoid chewing with the area where the surgery was performed. The possible administration of pre- and post-operative medicines remains at the discretion of the doctor. The stitches are to be removed after one week. Proper maintenance by the patient is recommended, including regular home oral hygiene and periodic check-ups in connection with professional hygiene sessions, which will extend the life of the device.

**SIDE EFFECTS:** After dental implant surgery, the most common side effects are pain, swelling, speech problems, and gingivitis. Other undesirable effects may result from pre-existing conditions. The dentist should inform the patient of possible common and specific undesirable effects in the case of pre-existing pathologies. Knowledge and identification of specific undesirable effects are a part of a dentist's study programme. The risks of implant surgery include but are not limited to: perforation of the labial or lingual plate, bone fractures, implant fractures, superstructure fractures, aesthetic problems, inadvertent perforation of bone and/or soft tissue (inferior alveolar canal, maxillary sinus, lingual artery, greater palatine foramen), nerve injuries, impairment of the natural dentition, failure or inadequate osseointegration, loss of bone or soft tissue volume due to concomitant events (local infection, peri-implantitis). Pathophysiological problems already diagnosed in the patient may increase the risks described above.

**DISPOSAL:** If it is necessary to dispose of medical devices, particularly when contaminated with blood or fluids, the appropriate hospital waste containers must be used. Use hospital grade gloves to avoid direct contact with the device in question.

**WARRANTY:** The C-Tech Implant Company offers a warranty on implants. In case of implant failure, the dentist should send the correct forms to [www.c-tech-implant.com/en/warranty](http://www.c-tech-implant.com/en/warranty) and to the dentist's distributor. Replacement implant(s) will be shipped to the dentist free of charge.

**INFORMATION FOR CUSTOMERS:** No person is authorised to provide information deviating from the information provided in this instruction sheet. For further information and technical documents, please consult our website at [www.c-tech-implant.com](http://www.c-tech-implant.com).

The Summary of Safety and Clinical Performance (SSCP) is available on EUDAMED at the following website:

<https://ec.europa.eu/tools/eudamed>. It is connected to the following BASIC UDI-DI:

- 805296102BIPHASICIMPLANWP (biphasic implants);

If any serious or suspected incidents occur as a result of the use of our medical devices, they should be reported immediately to the Company and to the national competent authority.

**"These instructions for use have been translated and are available at the following link: [www.c-tech-implant.com](http://www.c-tech-implant.com).**

**Alternatively, scan the following QRcode"**


In the event of a sudden malfunction of the company website due to force majeure, scan the QRcode to obtain the latest available revision. Alternatively, request the revision of the desired IFU at the following email address [quality@c-tech-implant.com](mailto:quality@c-tech-implant.com). The response will be sent as soon as possible. The paper version of the IFU in your language and which is among those available, can be requested from C-Tech Implant at the email address [quality@c-tech-implant.com](mailto:quality@c-tech-implant.com) at any time and will be sent immediately, compatibly with the working hours, and/or within a period not exceeding 7 working days. Alternatively, it can be requested at the time of the order and shipped at the same time, but not inside the packaging. This service is at no additional cost.

**DIRECT LINK TO INFORMATION MATERIALS:** [www.c-tech-implant.com](http://www.c-tech-implant.com)

**Chart 1**

Symbol	Description
	CE marking of a Medical Device with intervention of the Notified Body TÜV Rheinland
	Device Code
	Batch number
	Manufacturer
	Production Date
	Disposable Device
	Do not re-sterilise
	Expiry Date
	Warning
	Read the package leaflet carefully or consult the electronic IFU here: <a href="http://www.c-tech-implant.com">www.c-tech-implant.com</a>
	Do not use if packaging is damaged or opened
	Medical device sterilized using irradiation. Single sterile barrier system with protective packaging outside
	Keep away from direct sunlight
	Keep in a dry place
	Medical device
	Unique Device Identifier

**Chart 2**

Symbol	Description
	Patient Identification
	Information website for patient
	Nome of the doctor or the healthcare centre
	Date of the surgery
	Medical Device

**UDI Code and Data Matrix**

(01): UDI-DI, identifies the manufacturer and the specific version or model of device. (10): Batch number. (17): Expiry date. (11): Date of production.

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
Chart 1

Symbol	Description
	CE marking of a Medical Device with intervention of the Notified Body TÜV Rheinland
	Device Code
	Batch number
	Manufacturer
	Production Date
	Disposable Device
	Do not re-sterilise
	Expiry Date
	Warning
	Read the package leaflet carefully or consult the electronic IFU here: <a href="http://www.c-tech-implant.com">www.c-tech-implant.com</a>
	Do not use if packaging is damaged or opened
	Dispositivo médico esterilizado por irradiação. Sistema de barreira estéril único com embalagem protetora externa
	Keep away from direct sunlight
	Keep in a dry place
	Medical device
	Unique Device Identifier

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