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Product identification and intended use

Dental implants are medical devices made for use exclusively in the oral cavity, in patients suffering from total or partial mandibular and/or maxillary edentulism (without natural teeth). The devices are inserted into the bone to act as a substitute for the natural tooth root that will provide an anchor for the permanent dental crown.

They are implantable devices that come into contact with bone and soft tissue in the mouth. They have a long service life if the patient health conditions are ideal, and all the necessary maintenance interventions are carried out.

The external structure of an implant has a conical shape with a thread that allows it to make its way into harder bone, maintain its stability in softer bone and not damage the more sensitive areas of the mouth such as the canal nerve or Schneider's membrane (maxillary sinus). Implant threads are therefore different depending on the indications for use and the biological structures that patients present. The threads and apices can be more aggressive when they have to be inserted into a harder bone. On the other hand, they may be less aggressive if they have to be inserted into areas of softer bone or near more delicate areas (nerves, membranes).

All mini-implants are equipped with an internal connection system that allows the devices intended for prosthetic rehabilitation to be solidly anchored to the implants, thus allowing the restoration of missing elements. The prosthetic components of the SD and MB lines are compatible with each

other and can be used for all implant lengths according to the planned prosthetic rehabilitation. More details on the diameters and lengths available can be found under "Supply, Diameters and Lengths." Prosthetic devices may be available for both lines or only for one of them. It would be a good idea to consult the product catalogues to check availability for each line. Original components and products from C-Tech Implant must be used for prosthetic rehabilitations of the SD and MB lines. C-Tech disclaims any liability for the use of non-original instrumentation. It is recommended that only appropriate instruments supplied by C-Tech Implant, which are available individually or in special surgical organizers, be used for the surgical insertion procedure. C-Tech Implant disclaims any liability for the use of non-original equipment. For more details on the indications for use of each line, please see the section 'Indications and Contraindications'.

SD and MB implants are manufactured from grade 5 titanium and are subjected to a sandblasting surface treatment. The types of titanium used for production comply with harmonised standards.

IDENTIFICATION OF THE MANUFACTURER: The manufacturer of the SD and MB implant lines is:

C-TECH IMPLANT SRL, Via Cesare Battisti 2, 40123, Bologna, Tel. +39 051.6661817, email: info@c-tech-implant.com, website: www.c-tech-implant.com.

SUPPLY, DIAMETERS AND LENGTHS

Dental implants are supplied in a plastic vial inside a heat-sealed, sterile blister pack placed inside a cardboard box (secondary packaging). One box contains 1 unit. The blister pack and outer carton are duly labelled and indicate the code, expiry date, batch number and UDI code. The labels on the blister pack are removable and must be affixed by the medical personnel to be affixed to the patient implant card, which is also provided inside the implant case to then be delivered to the patient. They should also be attached to the patient's medical record. The Instructions for Use are also supplied with the implant and the patient implant card. The implant card must be filled in in all its fields according to what is requested and based on the information held by the medical staff. Please check the "Table 1" for a better understanding of the symbols on the medical device labels and the "Table 2" for a better understanding of the symbols on the patient implant card.

Implants are available in the following diameters and lengths:

LINE	LENGTHS AND DIAMETERS	
SD	O-Ball	Square head
	<u>Fine thread with collar</u> Ø1,8 L.10,11.5,13,15,18 mm Ø2,1 L.10,11.5,13,15,18 mm Ø2,5 L.10,13.15 mm <u>Aggressive thread with collar</u> Ø2,4 L.10,13.15,18 mm <u>Fine thread without collar</u> Ø1,8 L.10,13.15 mm Ø2,1 L.10,13.15 mm <u>Aggressive thread without collar</u> Ø2,4 L.10,13.15 mm	<u>Fine thread with collar</u> Ø1,8 L.10,13.15 mm Ø2,1 L.10,13.15 mm Ø2,5 L.10,13.15 mm <u>Aggressive thread with collar</u> Ø2,4 L.10,13.15 mm
MB	Mandibular Ø3.0 L. 9,11,13,15 mm	Mandibular Ø3.0 L.9,11,13,15 mm
	Maxillar Ø3.0 L. 9,11,13,15 mm	Maxillar Ø3.0 L. 9,11,13,15 mm

STORAGE AND USE: Implants should be stored in a cool, dry place at room temperature, away from heat, sunlight and water. Sterile products must be used no later than the expiry date stated on the labels. Do not use devices from open or damaged packages or blisters whose sterility can no longer be ensured.

WARNINGS: •C-Tech Implant disclaims all liability for damages caused to the patient by improper use or use other than as described in this instruction sheet or use by unqualified and/or untrained

personnel. • The law allows the surgical use of dental implants only and exclusively by qualified surgeons and dentists. Use by other persons is absolutely forbidden. • The manufacturer disclaims any liability for the use of unmarked medical devices not supplied by C-Tech Implant. • The surgical techniques required to place dental implants are highly specialised and complex procedures. Professionals will have to attend and complete courses to be ready for implantology techniques. • C-Tech Implant medical devices are provided with information materials and the company organises training courses worldwide. Participation in a specialised course is therefore strongly recommended. • The use of improper techniques can lead to patient injury, implant failure, bone loss and post-operative complications. • Dental implants are supplied in unopened, sterile packages and are ready for use. These are disposable and non-reusable devices. • Do not use a device in packaging that has been tampered with. • Do not re-use or re-sterilise the device. • Do not use the device beyond the expiry date. • Due to potential effects of anaesthesia do not use on pregnant women. • Use the devices on adults or on people who have reached the age of majority. • The device has not been evaluated in paediatric/adolescent patients and is not recommended for use on children. Routine treatment is not recommended until the end of the juvenile growth of the maxillary/mandibular bone has been properly documented. • Inform the healthcare professional of the presence of dental implants in case it is necessary to undergo magnetic fields, diagnostic investigations, therapeutic treatments or other procedures. In most cases, titanium dental implants do not interfere with magnetic fields and there are no known contraindications of any kind for therapeutic treatments or diagnostic investigations.

DIRECTIONS FOR USE:

SD Line: SD implants are ideal for those clinical situations where availability of the bone is reduced, whether vertically or horizontally. They can be used to stabilise removable prostheses by providing stability and attachment. They can be used to replace single missing tooth elements where implant space is limited or to rehabilitate just a few missing natural teeth. Bone quality permitting, they can be subjected to immediate loading. Otherwise, it will be necessary to wait for a period of time for integration. Four is the minimum number of mini-implants required to stabilise a lower arch prostheses. Six implants are recommended for upper arch prostheses. **Do not use:** If the conditions for dental implant surgery are not excellent and do not comply with the conditions outlined above. After a careful review of the patient's case, only the doctor can assess whether or not to use this line.

MB Line: MB implants are ideal for those clinical situations where availability of the bone is reduced, whether vertically or horizontally. They can be used to stabilise removable prostheses by providing stability and attachment. They can be used to replace single missing tooth elements where implant space is limited or to rehabilitate just a few missing natural teeth. Bone quality permitting, they can be subjected to immediate loading. Otherwise, it will be necessary to wait for a period of time for integration. Four is the minimum number of mini-implants required to stabilise a lower arch prostheses. Six implants are recommended for upper arch prostheses. **Do not use:** If the conditions for dental implant surgery are not excellent and do not comply with the conditions outlined above. After a careful review of the patient's case, only the doctor can assess whether or not to use this line.

The information provided by C-Tech Implant about each implant line should help the clinician in choosing the most appropriate device to use for each case. It is, however, the user's responsibility to make the final choice of the implant to be used once they have assessed its suitability for the case and based on their own experience and medical expertise.

CONTRAINDICATIONS: Established allergies to Ti-6AL-4V titanium alloy, recent myocardial infarction, valvular prostheses, severe renal disorders, severe kidney disorders, osteomalacia, treatment-resistant diabetes, recent high-dose radiotherapy, severe chronic alcoholism, drug use, terminal illnesses and tumours. **General:** Uncooperative patient. Intraoral contraindications such as: Anatomical impediments, insufficient bone tissue or poor bone quality. **Related:** Chemotherapy, mild renal disorders, mild liver disorders, endocrine disorders, psychological disorders or psychosis, lack of understanding and motivation, AIDS, seropositive cases, prolonged use of corticosteroids, calcium-phosphorus metabolism disorders, erythropoietic disorders.

WARNING! The success of an implant cannot be 100% guaranteed if the indications for use of the

product are not followed and the appropriate surgical/management procedures are not put in place. Close collaboration among the surgeon, the prosthetist and the dental laboratory is essential for successful implant treatment. When using a new device or treatment method for the first time, it is advisable to work alongside a colleague experienced in the use of the new device/treatment method to try to avoid possible complications. C-Tech Implant will gladly put you in contact with a product manager, its network of consultants and key opinion leaders.

PATIENT SELECTION: The medical devices covered by this leaflet are to be used on patients who have been diagnosed with total and/or partial mandibular and/or maxillary edentulism. Use the devices on adults or on people who have reached the age of majority. The device has not been evaluated in paediatric/adolescent patients and is not recommended for use on children. Routine treatment is not recommended until the end of the juvenile growth of the maxillary/mandibular bone has been properly documented. Ethnicity has no impact on the use of these medical devices. Not all edentulous patients are amenable to implant treatment, which is why it is necessary for medical personnel to subject the patient to an appropriate and careful initial selection process and then a precise diagnosis. Patient assessment prior to an implantology surgical procedure is extremely important. This includes assessment of the patient's general health, oral hygiene status, dental and periodontal condition and their anatomical acceptability. Taking the patient's medical and clinical history is essential. Allergies to titanium are rare but possible. It is advisable to check in advance with patients not have allergies of this type. The manufacturer provides the technical data sheets of the materials upon request for a deeper investigation. Particular attention should be paid to patients with local or systemic factors that could interfere with or slow down the healing processes of bone, soft tissue and osseointegration (e.g. smoking, poor oral hygiene, uncontrolled diabetes, orofacial radiotherapy, steroid therapy, infections in the surrounding bone). Special attention should be paid to patients undergoing bisphosphonate therapy. Patients with hyperfunction (bruxism) may be at greater risk of implant failure or fracture. In general, for patient selection and prosthetic design, the dentist must consider all the biological and physical conditions of the patient. The elements for this analysis are acquired during the doctor's studies and career path.

PRE-OPERATIVE PLANNING AND TREATMENT: Preparation for the surgical procedure includes: - General medical, psychological and dental history, general medical examination, clinical (complete blood count) and radiological examinations, CT scan and, if necessary, consultation with the patient's primary care doctor, - Patient information (indications, contraindications, clinical picture, possible undesirable effects, need for periodic follow-ups) – Dental hygiene plan, with any required periodontal procedures. -Evaluation of the procedure and surgical technique based on the treating dentist's studies, knowledge and experience. Since surgical techniques are highly specialised and complex procedures, it is highly recommended that practitioners complete a training course to learn established placement techniques in oral implantology. Failure to seek proper training or to follow the instructions for use (IFU) can affect the clinical success of the implant, such as failure, bone loss and postoperative complications, which could lead to pain for the patient and put his or her health at serious risk. Surgical placement and subsequent reconstruction with dental implants require complex procedures that may involve risks. It is the doctor's responsibility to inform the patient of any risks and alternatives (including a 'no treatment' option) before performing any procedure alternatives.

DIRECTIONS FOR USE AND SURGICAL PROCEDURE: Surgical techniques for implants are generally taught in university undergraduate programmes in dentistry. Nevertheless, the recommendations in this section should be considered. Users are encouraged to attend professional training courses and/or to read the information materials provided by the company and available on the website at www.c-tech-implant.com. At the bottom of this IFU there is a QR code that links to the page on the information material site. WARNING! The surgical procedure for inserting implants into the oral cavity must take place in a field that is as sterile as possible. Medical personnel are to be equipped with all necessary personal protective equipment (e.g., gloves, mask, gown, etc.) during surgery and whenever handling any medical devices. To proceed with the surgical procedure, medical personnel should be equipped with the surgical instrumentation sold by C-Tech Implant, which must be compatible with the implant line and the surgical technique to be used. The surgical instruments supplied by the manufacturer must be cleaned, disinfected and sterilised before first

use and after each re-use in accordance with the IFU sheet accompanying those instruments. For proper use of the surgical instrumentation and the surgical protocol to be put into practice, users are encouraged to attend professional training courses and/or to read the information materials provided by the company and available on the website at www.c-tech-implant.com. At the bottom of this IFU there is a QR code that links to the page on the information material site. The choice of the surgical procedure to be performed is to be assessed and decided upon by the doctor, in collaboration with the dental technician, after having carefully taken the patient's medical history and only once the prosthetic project has been properly planned and organised.

There are two surgical procedures that can be performed with C-Tech devices: traditional surgery and guided surgery. Following the evaluation of the patient and the corresponding radiographic overview, the type and number of implants to be placed will be determined. A minimum of four implants is recommended for mandibular cases and a minimum of six for the maxilla. The transparencies supplied with the surgical kit should be used to determine the length of the implant. A minimum distance of 5 mm must be maintained between implants to allow for the matrices to be slotted. In mandibular cases, implants should be placed starting from at least 5 mm from the anterior portion of the mental foramen. Mini-implants are self-tapping implants, which means that for insertion it is not necessary to drill a hole that is the same diameter as the implant itself. It is sufficient to under-prepare the implant site in a manner appropriate to the bone type.

Traditional surgery: The speed selected by the dentist when preparing the implant site depends on the preparation technique chosen. The wide range available allows the dentist to choose the most appropriate technique whilst remaining within values that do not lead to bone tissue damage. For additional details on the surgical procedure, please attend the courses organised by the company and read the materials available on the website.

SD: The drill selected for the initial perforation (\varnothing 1.1 D1/D2 bone or \varnothing 1.3 D2/D3 bone), equipped with external sterile irrigation, is placed on the selected sites and is lightly moved up and down until the cortical plane is perforated (400-800 RPM). In the case of hard bone, proceed with a \varnothing 1.5 mm drill (400-800 RPM). Unless the crestal bone is defective or too thin, there will be no need to enlarge the flap or incision. In the case of medium bone hardness, the drilling depth should be about 1/3 of the implant length. In the case of compact bone, the depth should be extended to about 3/4 of the length of the implant.

MB: Remove a small portion of gingival tissue with a tissue punch and proceed to create an pilot hole for drilling with a locator drill. For the cutting speed, follow what is stated in this leaflet and in the surgical protocol. For the preparation of the implant site, it is necessary to start by creating a hole with a locator drill (400-800 RPM). Then proceed using the pilot drill (400- 800 RPM), with the desired stop length, to drill a hole of the length needed for the implant selected. In the case of hard bone, make one pass with the \varnothing 2.6mm drill (400-800 RPM). In the case of medium bone hardness, the drilling depth should be about 1/3 of the implant length. In the case of compact bone, the depth should be extended to about 3/4 of the length of the implant.

For additional details on the surgical procedure, please attend the courses organised by the company and read the materials available on the website.

Guided surgery: The speed selected by the dentist when preparing the implant site depends on the preparation technique chosen. The wide range available allows the dentist to choose the most appropriate technique whilst remaining within values that do not lead to bone tissue damage. For additional details on the surgical procedure, please attend the courses organised by the company and read the materials available on the website.

To successfully use guided surgery, the dentist must have good previous experience in traditional surgery and also very good knowledge of anatomy. For proper patient selection for this type of surgical procedure, follow the instructions given in the chapter on 'patient selection' and apply what it's learned during your studies. Remove the soft tissue portion using a tissue punch (100-800RPM). Then proceed by drilling a pilot hole for the drills with a locator drill (400-800 RPM). Then the dentist should proceed to drill sequentially (400-800 RPM), using as many drills at as many diameters as necessary, until reaching the diameter and drill length that correspond to the chosen implant length.

POSITIONING THE IMPLANTS:

The implant is removed from the sterile vial and the implant tip is inserted into the implant site, using the cap as the initial wrench. The implant is pushed into the bone by exerting downward pressure and turning it clockwise. Once the resistance of the bone becomes too strong to proceed further, release the plug from the implant and continue with the butterfly driver to apply more force. In some cases the implant can be inserted completely with this last driver. However, if the insertion presents too much resistance, it will be necessary to use the ratchet wrench for final positioning.

Screw the implant down into the bone with the contra-angle torque wrench (at 35 Ncm up to a maximum of 45 Ncm), then as required, use the ratchet torque wrench (at 35 Ncm up to a maximum of 45 Ncm), turning the implant clockwise until it has been completely positioned in the bone. If, in order to achieve complete insertion, greater force is required, the implant must be carefully unscrewed and removed. Then use the pilot drill again to make the implant site deeper.

POST-OPERATIVE COURSE: The patient should be supervised and instructed to follow a post-operative course appropriate to the surgery undergone. This includes the application of cold compresses, and no physical exertion and/or sauna for the first 24 hours. Alcohol, nicotine, tea or coffee should be limited for the first few days, and the patient should avoid chewing with the area where the surgery was performed. Proper maintenance by the patient is recommended, including regular home oral hygiene and periodic check-ups in connection with professional hygiene sessions, which will extend the life of the device.

SIDE EFFECTS: After dental implant surgery, the most common side effects are pain, swelling, speech problems, and gingivitis. Other undesirable effects may result from pre-existing conditions. The dentist should inform the patient of possible common and specific undesirable effects in the case of pre-existing pathologies. Knowledge and identification of specific undesirable effects are a part of a dentist's study programme. The risks of implant surgery include but are not limited to: perforation of the labial or lingual plate, bone fractures, implant fractures, superstructure fractures, aesthetic problems, inadvertent perforation of bone and/or soft tissue (inferior alveolar canal, maxillary sinus, lingual artery, greater palatine foramen), nerve injuries, impairment of the natural dentition, failure or inadequate osseointegration, loss of bone or soft tissue volume due to concomitant events (local infection, peri-implantitis). Pathophysiological problems already diagnosed in the patient may increase the risks described above.

DISPOSAL: If it is necessary to dispose of medical devices, particularly when contaminated with blood or fluids, the appropriate hospital waste containers must be used, using hospital grade gloves to avoid direct contact with the device in question.

WARRANTY: The C-Tech Implant Company offers a warranty on implants. In case of implant failure, the dentist should send the correct forms to www.c-tech-implant.com/en/warranty and to the dentist's distributor. Replacement implant(s) will be shipped to the dentist free of charge.

INFORMATION FOR CUSTOMERS: No person is authorised to provide information deviating from the information provided in this instruction sheet. For further information and technical documents, please consult our website www.c-tech-implant.com in the 'download' page.

The Summary of Safety and Clinical Performance (SSCP) is available on EUDAMED at the following website: <https://ec.europa.eu/tools/eudamed>. It is connected to the following BASIC UDI-DI'S:

- 805228249MONOPHASIMPLANB6 (monophasic implants);

If any serious or suspected incidents occur as a result of the use of our medical devices, they should be reported immediately to the Company and to the national competent authority.

"These instructions for use have been translated and are available at the following link: www.c-tech-implant.com . Alternatively, scan the following QRcode"

In the event of a sudden malfunction of the company website due to force majeure, scan the QRcode to obtain the latest available revision. Alternatively, request the revision of the desired IFU at the following email address quality@c-tech-implant.com . The response will be sent as soon as possible. The paper version of the IFU in your language and which is among those available, can be requested from C-Tech Implant at the email address quality@c-tech-implant.com at any time and will be sent immediately, compatibly with the working hours, and/or within a period not

exceeding 7 working days. Alternatively, it can be requested at the time of the order and shipped at the same time, but not inside the packaging. This service is at no additional cost.

DIRECT LINK TO INFORMATION MATERIALS: www.c-tech-implant.com

Chart 1

Symbol	Description
	CE marking of a Medical Device with intervention of the Notified Body TÜV Rheinland
	Device Code
	Batch number
	Manufacturer
	Production Date
	Disposable Device
	Do not re-sterilise
	Expiry Date
	Attention
	Read the package leaflet carefully or consult the electronic IFU here: www.c-tech-implant.com
	Do not use if packaging is damaged or opened
	Medical device sterilized using irradiation. Single sterile barrier system with protective packaging outside
	Keep away from direct sunlight
	Keep in a dry place
	Medical device

UDI	Unique Device Identifier
------------	--------------------------

Chart 2

Symbol	Description
	Patient Identification
	Information website for patient
	Nome of the doctor or the healthcare centre
	Date of the surgery
	Medical Device

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