EN These instructions for use have been translated and are available at the following link: [www.c-tech-implant.com](http://www.c-tech-implant.com). Otherwise, scan the following QRcode. Alternatively, request a review of your desired IFU at the following email address [quality@c-tech-implant.com](mailto:quality@c-tech-implant.com). For further information, read the submission methods in the "Information for Customers" paragraph.

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**(EN) DENTAL DRILLS AND COMPONENTS FOR DENTAL DRILLS**

**PRODUCT IDENTIFICATION AND INTENDED USE**

Dental drills and components are medical devices made for the exclusive use in the oral cavity, in patients suffering from total or partial edentulism (absence of teeth). They are involved in the dental implant site preparation process inside the oral cavity. Dental drills and components for dental drills are medical devices that are meant to be used exclusively with C-Tech Implant Dental implants. The use of non-original components limits the responsibility of C-Tech Implant Srl on the good outcome of the surgical procedure. The clinical benefit derived from the use of these devices is linked to the one of the dental implants and is aesthetic and functional improvement in relation to the patient's clinical situation (edentulism) resulting in improved quality of life of the patient.

**IDENTIFICATION OF THE MANUFACTURER:** The manufacturer of dental drills and components for dental drills is: **C-TECH IMPLANT SRL** Via Cesare Battisti 2, 40123, Bologna, Tel. +39 051.6661817, email: [info@c-tech-implant.com](mailto:info@c-tech-implant.com), website: [www.c-tech-implant.com](http://www.c-tech-implant.com)

**SUPPLY AND MODELS**

Dental Drills and components for dental drills are supplied in sealed polypropylene pouches in non-sterile conditions. Each pouch is supplied with the label indicating the medical device reference, the batch number and UDI code. Each medical device is also supplied with the instruction for use. The devices can be also supplied in boxes that are used as organizers for the devices during the surgical operation. For additional details on the devices models and availability of each implant line (EL,CP,ND, BL,GL,SD and MB), please refer to the product catalog of each line in this section of the website [www.c-tech-implant.com](http://www.c-tech-implant.com) or by scanning the QR code at the end of this IFU. For first use and reprocessing procedures please follow the instructions in paragraph "INSTRUCTIONS FOR CLEANING, DRYING AND STERILISATION".

**WARNINGS:** C-Tech Implant declines all responsibility for damage caused to the patient by improper use or different use from what is reported in this instruction sheet or for use by unqualified and / or untrained personnel. The law allows the use of medical devices only and exclusively to qualified surgeons and dentists. Use by other people is strictly prohibited. • C-Tech Implant has no responsibility for the use of unmarked medical devices not supplied by the manufacturer. • For the surgical procedure of each implant line refer to the specific surgical protocol with the identification of the correct instruments to be used. • C-Tech Implant medical devices are equipped with informative material (IFU and catalogs) and the company organizes training courses all over the world to give the necessary information regarding the use of its medical devices. Participation in a specialization course is therefore strongly recommended. • The use of improper techniques can cause patient harm, implant failure, bone loss and post-operative complications. • Do not reuse or re-sterilize disposable devices. • The maximum number of reuses of each device is 30 times. Any use and reuse of the surgical instrumentation without a prior careful and in-depth reading of this package leaflet or of the catalogues indicated below will affect the functional performance and the safety of the devices provided for by the manufacturer, exposing the patient to the risk of a failure of the product and/or to serious risks of contamination. • Do not use a device in packaging that has been tampered with. • Due to potential effects of anesthesia do not use on pregnant women. • Use the devices on adults or on people who have reached the age of majority. • The device has not been evaluated in paediatric/adolescent patients and is not recommended for use on children. Routine treatment is not recommended until the end of the juvenile growth of the maxillary/mandibular bone has been properly documented.

**STORAGE AND USE:** Store in a cool, dry place at room temperature, away from heat, sunlight and water. Devices that have undergone the reprocessing phase according to what is indicated in paragraph (INSTRUCTIONS FOR CLEANING, DRYING AND STERILISATION) must not be used after 6 months after the sterilization. Do not use devices from open or damaged packages or pouches whose sterile status can no longer be ensured.

For disposable devices, see the details in the paragraph "indications for use". C-Tech Implant declines all responsibility for the reuse of these devices on multiple patients.

**INDICATION FOR USE:**

**PLEASE NOTE: 1)** All the devices object of this instruction for use are supplied as not sterile devices. Before use and after each use,

they must be cleaned, disinfected and sterilized in accordance with the processes indicated in paragraph "INSTRUCTIONS FOR CLEANING, DRYING AND STERILISATION". The devices can be reused up to 30 times.

**2)** The speed selected by the dentist when preparing the implant site depends on the preparation technique chosen. The wide range available allows the dentist to choose the most appropriate technique whilst remaining within values that do not lead to bone tissue damage. For additional details on the surgical procedure, please attend the courses organized by the company and read the materials available on the website.

**DENTAL DRILLS**

- A) LOCATOR DRILL- Description:** They are used to create an incision on the cortical bone to facilitate the work of the following instrument. **Material and Packaging:** Manufactured in stainless steel, supplied in non-sterile packaging. The device is reusable. **RPM:** 100-750 RPM
- B) COUNTERBORE-Description:** They are used to reduce pressure around the cortical bone to facilitate the accommodation of the dental implant shoulder. **Material and Packaging:** Manufactured in stainless steel, supplied in non-sterile packaging. The device is reusable. **RPM:** 200-600 RPM
- C) PREPARATION DRILLS- Description:** They are used to create the hole of the implant of the desired diameter and length. **Material and Packaging:** Manufactured in stainless steel and titanium grade 5, supplied in non-sterile packaging. The device is reusable. One dental drill of SD and MB line is sold as disposable. **RPM:**100-750 RPM (traditional surgery); 400-800 RPM (SD drills and guided surgery)
- D) TISSUE PUNCH- Description:** It's used to remove a portion of the soft tissue in correspondence of the planned position of the implant. **Material and Packaging:** Manufactured in stainless steel, supplied in non-sterile packaging. The device is reusable. **RPM:** 100-800 RPM
- E) BONE PROFILER (FOR GUIDED SURGERY)- Description:** It's used to flatten the irregularities of the crestal bone **Material and Packaging:** Manufactured in stainless steel, supplied in non-sterile packaging. The device is reusable. **RPM:** 150-400 RPM
- F) BONE PROFILER (TRADITIONAL SURGERY)- Description:** It's used to remove bone around the coronal area of the implant to allow proper seating of Multi unit abutments. **Material and Packaging:** Manufactured in stainless steel, supplied in non-sterile packaging. The device is reusable. **RPM:** 150-400 RPM
- G) DRILL FOR PINS- Description:** It's used to create a hole for the accommodation of the lateral pins. **Material and Packaging:** Manufactured in stainless steel, supplied in non-sterile packaging. The device is reusable. **RPM:** 150-800 RPM
- H) BONE EXPANDER- Description:** It's used to compact bone inside the designed implant site. **Material and Packaging:** Manufactured in titanium grade 5, supplied in non-sterile packaging. The device is reusable. **RPM:** 15-40 RPM
- I) BONE TAPS- Description:** Used to thread the bone inside the designed implant site. **Material and Packaging:** Manufactured in stainless steel, supplied in non-sterile packaging. The device is reusable. **RPM:** 15-30 RPM

**COMPONENTS FOR DENTAL DRILLS**

- A) METAL DEPTH STOPPER- Description:** Connected to the drills allows the perforation of the implant site of the desired length. **Material and Packaging:** Manufactured in stainless steel, supplied in non-sterile packaging. The device is reusable.
- B) EXTENDER- Description:** connected to the surgical instrument allows to extend the length of the instruments to facilitate their use in specific narrow areas. **Material and Packaging:** Manufactured in stainless steel, supplied in non-sterile packaging. The device is reusable.
- C) PARALLEL PINS: Description:** Used to check if the implant sites are drilled in a parallel way. **Material and Packaging:** Manufactured in stainless steel, supplied in non-sterile packaging. The device is reusable.
- D) DEPTH GAUGE- Description:** used to check the depth of the drilled implant site. **Material and Packaging:** Manufactured in stainless steel, supplied in non-sterile packaging. The device is reusable.
- E) LATERAL PIN- Description:** used to hold the surgical guide into position. **Material and Packaging:** Manufactured in stainless steel, supplied in non-sterile packaging. The device is reusable.
- F) ALL ON FOUR GUIDE ARCH- Description:** used to give a reference of the angulation of the multi-unit abutment (17° and 30°) during the drill phase of the implant site. **Material and Packaging:** Manufactured in Titanium grade 1 and 5, supplied in non-sterile packaging. The device is reusable.
- G) PIN FOR BONE PROFILER: Description:** used to guide the insertion of the bone profiler inside the bone for drilling. **Material and Packaging:** Manufactured in Titanium, supplied in non-sterile packaging. The device is reusable.
- H) BOX - Description:** The box allows you to collect and organize the medical devices inside. **Material and supply:** Manufactured from polymeric material and aluminium, supplied in non-sterile packaging but sterilized by the user. The devices are reusable. For the cleaning and sterilization procedure before use, see the chapter "CLEANING, DRYING AND STERILIZATION INSTRUCTIONS".

**CONTRAINDICATIONS:** The contraindication are strictly related to the contraindications of the dental implant system. **General:** Uncooperative patient. Intraoral contraindications such as: Anatomical impediments, insufficient bone tissue or poor bone quality. **Related:** Chemotherapy, mild renal disorders, mild liver disorders, endocrine disorders, psychological disorders or psychosis, lack of understanding and motivation, AIDS, seropositive cases, prolonged use of corticosteroids, calcium-phosphorus metabolism disorders, erythropoietic disorders.

Additional contraindications are known allergies to titanium grade 1 and 5, and stainless steel.

**PATIENT SELECTION:** Patient selection phase is done during the initial examination of the patient for dental implant suitability.

The medical devices covered by this leaflet are to be used on patients who have been diagnosed with total and/or partial mandibular and/or maxillary edentulism. The device has not been evaluated in paediatric/adolescent patients and is not recommended for use on children. Routine treatment is not recommended until the end of the juvenile growth of the maxillary/mandibular bone has been properly documented. Ethnicity has no impact on the use of these medical devices.

All the raw materials are widely used in the medical field and biocompatible, allergies are rare but possible. It is advisable to check in advance with patients not to have allergies of this type. The manufacturer provides the technical data sheets of the materials upon request for a deeper investigation.

**PRE-OPERATIVE PLANNING AND TREATMENT:** The suitability of a patient, affected by total or partial edentulism, for implant surgery is derived from a profound analysis of the patient and his clinical history. This analysis is carried out by the surgeon and the dental technician on the basis of the skills acquired during their studies and experience. Once the suitability of the patient for dental implantology has been defined, it is always the task of the surgeon and the dental technician to carry out the surgical and prosthetic planning based on the chosen approach.

**POST-OPERATIVE RECOVERY:** The patient must be informed about the post-operative recovery following surgery.

**SIDE EFFECTS:** After dental implant surgery, the most common side effects are pain, swelling, speech problems, and gingivitis. Other undesirable effects may result from pre-existing conditions. The dentist should inform the patient of possible common and specific undesirable effects in the case of pre-existing pathologies. Knowledge and identification of specific undesirable effects are a part of a dentist's study program. The risks of implant surgery include but are not limited to: perforation of the labial or lingual plate, bone fractures, implant fractures, superstructure fractures, aesthetic problems, inadvertent perforation of bone and/or soft tissue (inferior alveolar canal, maxillary sinus, lingual artery, greater palatine foramen), nerve injuries, impairment of the natural dentition, failure or inadequate osseointegration, loss of bone or soft tissue volume due to concomitant events (local infection, peri-implantitis). Pathophysiological problems already diagnosed in the patient may increase the risks described above.

**DISPOSAL:** If it is necessary to dispose of medical devices, particularly when contaminated with blood or fluids, the appropriate hospital waste containers must be used. Use hospital grade gloves to avoid direct contact with the device in question.

#### **INSTRUCTIONS FOR CLEANING, DRYING AND STERILISATION**

It is recommended to strictly follow the instructions below for a correct and effective cleaning, disinfection and sterilization procedure.

**PLEASE NOTE:** The medical devices supplied by C-Tech Srl as "non-sterile", must be cleaned and sterilised before their first use. The cleaning and sterilisation must be performed according to the methods indicated in the paragraphs below. Reusable devices, before the first use and for all the subsequent uses (max.30 times), must always be subjected to all the cleaning, disinfection, drying and sterilisation steps described below. It is recommended that cleaning operations be carried out as soon as the surgical procedure is completed to minimize the risk of contaminants drying on the devices, making their removal more difficult.

**WARNINGS** • *Strictly follow the steps described below and carry them out sequentially, without skipping any step.* • *Always use, at all stages, personal protective equipment (such as gloves, mask, goggles, etc.).* • *The devices must not be placed in contact with other objects during the washing phases.* • *Devices that are contained in boxes must be taken out of them and washed separately. The latter must be reinserted in the box, after any reassembly, only after completing all the cleaning phases.* • *Do not use wire brushes or steel wool for cleaning. Residues adhering to instruments, drying up, can cause corrosion.* • *Never use chlorine-based disinfectants.* • *Chemical sterilization is not recommended as this procedure can damage the plastic surface of the surgical boxes and can oxidize the devices.* • *Do not use dry sterilization as the high temperature (about 180 ° C) could melt the plastic boxes.* • *Do not sterilize in glass ball sterilizers (temperature above 300 ° C).* • *Subject the autoclave, ultrasonic bath and the automated disinfectant to regular maintenance and checks as required by the user manual provided by the manufacturer of these machines.*

#### **CLEANING STEPS- MANUAL OPTION**

**1. PRE-TREATMENT:** Disassemble each part of the multiple instruments (when applicable). Remove the coarsest impurities from the instruments by rinsing with running water for at least 2 minutes. Then remove small residues of blood, saliva, tissue or bone from the instruments using a soft brush with running water, for at least 2 minutes. Do not use brushes with metal or steel wool bristles. The residues adhering to the instruments may cause corrosion as they dry. Such operation must be performed no more than 2 hours after use.

**2. CLEANSING:** Immediately after rinsing with running water (pre-treatment), immerse disassembled and dirty instruments in a basin with detergent solution for at least one hour, in order to prevent blood residues from coagulating. The instruments must be completely immersed in the solution and the basin must be stirred several times to make the process easier. It is recommended to use a solution based on peracetic acid, such as "SEKUSEPT™ PULVER CLASSIC" with a concentration of 2%. In particular, thoroughly clean the inner side of the cavities of the instruments.

**3. DISINFECTION:** For a thorough cleaning, place the instruments in the ultrasound machine and leave them immersed for at least 30 minutes and up to an hour. To prevent the corrosion of surgical instruments, such liquids must be suitable for the disinfection of metal and plastic instruments. It is recommended to use a disinfectant based on peracetic acid, such as "SEKUSEPT™ PULVER CLASSIC" with a concentration of 4%. Never use chlorine-based disinfectant. Wash container and instruments separately. Complete the preparation step of steam sterilisation with a thorough and careful rinsing, with sterile water, and then dry as described in the paragraph below.

**4.DRYING STEP:** Dry the clean instruments in a special hot-air oven at 60 °C for 2 hours. It is recommended to carefully follow the instructions for use provided by the manufacturer of the oven. Make sure that instruments are completely dry, since any residue of steam might start corrosive processes which damage the instrument over time.

#### **AUTOMATED CLEANING & DISINFECTION OPTION**

To carry out automated cleaning, disinfection and drying, follow the instructions provided by the manufacturer of the washer disinfectant in use. C-Tech Implant validated the cycle with a washer disinfectant model AT-OS AWD655-10 (DIN 10) that is programmable and equipped with data integrity validated SW. The process is described below.

•Place the devices and any surgical box in the appropriate basket suitable for mechanical washing, close it and load it into the washer disinfectant. •It is recommended to avoid contact of the devices with each other during power washing. •Start the appropriate washing, disinfection and drying programme for the devices. •It requires a pre-rinse with demineralised water for 3 minutes. •Cleaning should be carried out for 5 minutes at a temperature of 55°C, followed by rinsing for 2 minutes at a temperature of 40°C. •We recommend using the concentrated liquid alkaline detergent Neodisher MediClean Forte (0.5%) and the alkaline residue neutraliser Neodisher Z (0.1%). •Final rinse with demineralised water for 2 minutes. •The cycle continues with automated disinfection at 90°C for 5 min. •Take the devices from the washer disinfectant at the end of the drying program and proceed to the sterilisation phase.

**STERILISATION STEPS 1) PACKAGING BEFORE STERILISATION** -The disinfected, clean and dry instruments must be tidily arranged in the surgical box, which must be sealed in special sterilisable envelopes validated according to ISO 11607-1 and according to what is indicated in the instruction manual supplied by the manufacturer. Alternatively, each instrument can be placed inside paper envelopes (60 g/m<sup>2</sup>) and plastic film (PET/PP). The manufacturer must indicate that the paper envelopes (60 g/m<sup>2</sup>) and plastic film (PET/PP) can be steam-sterilised. The paper envelopes and plastic film must be sealed onto all four sides. C-Tech recommends using the same sealing machine model and the same sealing parameters used during the validation phase: LUX SEAL M9025/4 sealing machine; sealing time

of 5 seconds; sealing temperature: 150 °C. **2) STERILISATION** For a long service life of the surgical box, sterilisation must be performed in autoclave with steam at a temperature of 121 °C for at least 15 minutes. After sterilization let dry for 20 minutes. It is recommended to validate your autoclave with full and worst-case load in accordance with standard ISO 17665-1. Store the box or the devices individually packaged in a dry environment and inside the relevant envelope for a maximum of 6 months from the sterilisation date. C-Tech Implant recommends using an autoclave model equivalent to the one used during the validation phase: "SERENA 18" model (Manufacturer: Reverberi Srl).

**INFORMATION FOR CUSTOMERS:** No person is authorised to provide information deviating from the information provided in this instruction sheet. For further information and technical documents, please consult our website at [www.c-tech-implant.com](http://www.c-tech-implant.com). If any serious or suspected incidents occur as a result of the use of our medical devices, they should be reported immediately to the Company and to the national competent authority.

DIRECT LINK TO INFORMATION MATERIALS: [www.c-tech-implant.com](http://www.c-tech-implant.com)

**If any serious or suspected incidents occur as a result of the use of our medical devices, they should be reported immediately to the Company and to the national competent authority.**

**"These instructions for use have been translated and are available at the following link: [www.c-tech-implant.com](http://www.c-tech-implant.com).**

**Alternatively, scan the following QRcode".**

In the event of a sudden malfunction of the company website due to force majeure, scan the QRcode to obtain the latest available revision. Alternatively, request the revision of the desired IFU at the following email address [quality@c-tech-implant.com](mailto:quality@c-tech-implant.com). The response will be sent as soon as possible. The paper version of the IFU in your language and which is among those available, can be requested from C-Tech Implant at the email address [quality@c-tech-implant.com](mailto:quality@c-tech-implant.com) at any time and will be sent immediately, compatibly with the working hours, and/or within a period not exceeding 7 working days. Alternatively, it can be requested at the time of the order and shipped at the same time, but not inside the packaging. This service is at no additional cost.

### SYMBOLS ON THE LABELS OF MEDICAL DEVICES

Table 1: Definition of the symbols on the labels.

Symbol	Description
	CE Marking of a Medical Device If applicable: the identification number of the notified body must follow this symbol.
	CE marking of a Medical Device with intervention of the Notified Body TÜV Rheinland
	Code of the device
	Lot Number
	Manufacturer
	Manufacturing date
	Non Sterile device
	Disposable device
	Warning
	Read the package leaflet carefully or consult the electronic IFU here: <a href="http://www.c-tech-implant.com">www.c-tech-implant.com</a>
	Do not use in case of damaged packaging
	Keep away from direct sunlight
	Keep in a dry place
	Medical device
	Unique Device Identifier

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