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COMPONENTS FOR DENTAL IMPLANTS

PRODUCT IDENTIFICATION AND INTENDED USE: Components for dental implants are medical devices made for the exclusive use in the oral cavity, in patients suffering from total or partial edentulism (absence of teeth), mandibular and / or maxillary where, to treat this clinical condition, one or more dental implants have been inserted. They are involved during the phases of planning and creation of the final prosthesis. The components for implants are medical devices that are meant to be used in combination with C-Tech Dental implants. The interconnections between dental implants and components or among components, is easily understandable thanks to the information in the product catalogs, the use of codes with references to the implant lines and the constant formation organized by the company.

The use of non-original components limits the responsibility of C-Tech Implant Srl on the good outcome of the surgical and prosthetic procedure. Prosthetic components must be inserted and placed inside the implants with the appropriate tools. It is recommended to use tools manufactured by C-Tech Implant.

Clinical benefits of the devices for the patient are: aesthetic and functional improvement in relation to the patient's clinical situation (edentulism) resulting in improved quality of life of the patient, a long service life if the patient health conditions are ideal, and all the necessary maintenance interventions are carried out. The clinical benefit of the components for dental implants is strictly related to the dental implants one. The components contribute to achieving a the specific benefit, in combination with the dental implants.

IDENTIFICATION OF THE MANUFACTURER: The manufacturer of components for dental implant is:

C-TECH IMPLANT SRL, Via Cesare Battisti 2, 40123, Bologna, Tel. +39 051.6661817, email: info@c-tech-implant.com

website: www.c-tech-implant.com

SUPPLY AND MODELS

Components for dental implants are supplied in sealed polypropylene pouches in sterile or non-sterile conditions depending on the type of medical device see details in paragraph "Indications for use". Each pouch is supplied with the label indicating the medical device reference, the batch number and UDI code. Each medical device is also supplied with the instruction for use. Some devices can be also supplied in boxes that are used as organizers for the devices during the surgical operation or the prosthetic restoration.

For additional details on the devices models and availability of each implant line (EL, CP, ND, BL, GL, SD and MB), please refer to the product catalog of each line in this section of the website www.c-tech-implant.com.

The implant card is supplied inside the dental implant packaging. For the components, the implant card must be filled in in the applicable field according to what is requested and based on the information held by the medical staff, also in case

of replacement of the component. Please check the "Table 1" for a better understanding of the symbols on the medical device labels and the "Table 2" for a better understanding of the symbols on the patient implant card.

For initial processing procedures please follow the instructions in paragraph "INSTRUCTIONS FOR CLEANING, DRYING AND STERILISATION".

WARNINGS: C-Tech Implant declines all responsibility for damage caused to the patient by improper use or different use from what is reported in this instruction sheet or for use by unqualified and / or untrained personnel. The law allows the use of medical devices only and exclusively to qualified surgeons and dentists. Use by other people is strictly prohibited.

• C-Tech Implant has no responsibility for the use of unmarked medical devices not supplied by the manufacturer. • Use only dental components related to the chosen implant line. •C-Tech Implant medical devices are equipped with informative material (IFU and catalogs) and the company organizes training courses all over the world to give the necessary information regarding the use of its medical devices. Participation in a specialization course is therefore strongly recommended. •The use of improper techniques can cause patient harm, implant failure, bone loss and post-operative complications. •Do not reuse or re-sterilize disposable devices. • Any use and reuse of the surgical instrumentation without a prior careful and in-depth reading of this package leaflet or of the catalogues indicated below will affect the functional performance and the safety of the devices provided for by the manufacturer, exposing the patient to the risk of a failure of the product and/or to serious risks of contamination. • Do not use a device in packaging that has been tampered with. •Due to potential effects of anaesthesia do not use on pregnant women. • Use the devices on adults or on people who have reached the age of majority. • The device has not been evaluated in paediatric/adolescent patients and is not recommended for use on children. Routine treatment is not recommended until the end of the juvenile growth of the maxillary/mandibular bone has been properly documented. • Inform the healthcare professional of the presence of dental implants and/or components in case it is necessary to undergo magnetic fields, diagnostic investigations, therapeutic treatments or other procedures. In most cases, titanium, cobalt chrome dental implants and/or components do not interfere with magnetic fields and there are no known contraindications of any kind for therapeutic treatments or diagnostic investigations. The cobalt chrome may interfere with final imaging resolution.

STORAGE AND USE: Store in a cool, dry place at room temperature, away from heat, sunlight and water. Sterile sold devices must be used no later than the expiry date stated on the labels. Devices that have undergone the reprocessing phase according to what is indicated in paragraph (INSTRUCTIONS FOR CLEANING, DRYING AND STERILISATION) must not be used after 6 months after the sterilization. Do not use devices from open or damaged packages or pouches whose sterile status can no longer be ensured.

Some components for dental implants are disposable devices, see the details in the paragraph "indications for use". C-Tech Implant declines all responsibility for the reuse of these devices on multiple patients.

INDICATIONS FOR USE

A) COVER SCREWS-Description: the cover screw is screwed into the implant once positioned in the implant site in order to allow its healing and avoiding possible inclusions of organic matter and bacteria inside the implant cavity.

Material and Packaging: Manufactured in grade 5 titanium, supplied in sterile packaging with the implant and also sold sterile in single pouches. The device is disposable. **Tightening Torque:** 10Ncm **Note:** Its use is always necessary unless the prosthetic protocol that the end user decides to implement provides for the use of other healing components or other prosthetic solutions.

B) HEALING ABUTMENTS: HEALING SCREWS-Description: the healing screw is screwed into the abutment to prevent the infiltration of material to allow healing and conditioning of soft tissues to ensure adequate housing of the chosen prosthesis. **Material and Supply:** Manufactured in titanium grade 5 and Peek, supplied in sterile packaging. The device is disposable. **Torque tightening:** 15Ncm **Note:** Its use is always necessary unless the prosthetic protocol that the end user decides to implement provides for the use of other healing components or other prosthetic solutions.

The lifetime of the device depends on the osseointegration of the dental implant which is evaluated by the clinician on each case.

HEALING CAP- Description: The healing cap is screwed onto the abutment to prevent the infiltration of material and allows healing and conditioning of the soft tissues waiting for the load of the final prosthesis. **Material and Supply:** Manufactured in grade 5 titanium, supplied in sterile packaging. The device is disposable. **Torque tightening:** 10Ncm **Note:** Its use is always necessary unless the prosthetic protocol that the end user decides to implement provides for the use of other healing components or other prosthetic solutions.

The lifetime of the device depends on the osseointegration of the dental implant which is evaluated by the clinician on each case.

C) TRANSFERS AND SCAN BODY-Description: the devices for taking impressions (traditional or digital) serve as reference points for the placement of devices in the mouth to transfer their position into an analog or digital reality, in order to reproduce models of the mouth on which dental technicians can work on to make the prosthesis. **Material and Supply:** Manufactured in titanium grade 5, Peek and Delrin in sterile packaging. The devices are disposable. **Torque tightening:** The Open and Closed tray Transfers must be screwed by hand, no indication of tightening torque is required. The scan caps must be tightened with screw at 25Ncm. The scanbody should be tightened with screw at 15Ncm. **Note:** Its use is always necessary and the choice between the various models available and the two solutions, analog or digital, depends on the prosthetic protocol that the end user decides to implement and from his own experience and medical skills.

D) ABUTMENTS-Description: The dental implant abutment is screwed to the dental implant and serves as the basis for the design and construction of the temporary or definitive prosthesis. **Material and Supply:** Manufactured in Titanium grade 5, Peek, Plexiglass, Chrome Cobalt. Supplied in non-sterile packaging but sterilizable by the user. The device is disposable. For before use cleaning and sterilization procedure see chapter "INSTRUCTIONS FOR CLEANING, DRYING AND STERILISATION". **Torque tightening:** For tightening indications see what is reported in the prosthetic screw section.

E) PROSTHETIC SCREWS- Description: The prosthetic screw is used to connect the prosthetic part to the implant, the prosthesis to the abutment or the impression taking devices to the implant or abutment. **Material and Supply:**

Manufactured in grade 5 titanium, supplied in non-sterile packaging but sterilized by the user. For before use cleaning and sterilization procedure see chapter "INSTRUCTIONS FOR CLEANING, DRYING AND STERILISATION". Screws for transfers are sold sterile. The devices are disposable. **Torque tightening:** 25 Ncm, if the prosthesis is to be fixed on implants that have a diameter equal to or greater than 3.8 mm; 20 Ncm, if the prosthesis is to be fixed on implants with a diameter of less than 3.8 mm and on OMNI abutments, 15Ncm for ND scanbodies; 15 Ncm in case of tightening of on "all-on-four" or "all-on-six" prosthesis, bars or bridges screwed on MUA/Flat abutments. Hexalobed screws: 20Ncm for Dual Channel abutments, 15Ncm for FLAT, MUA and OMNI abutments.

Class I Medical Devices

F) OBALL HOUSING-Description: The o-ball metal housing is placed inside the prosthesis and its meant to connect it to the o-ball implants or abutments. **Material and Supply:** Manufactured in grade 5 titanium (with o-ring or plastic cap) with different retentions, supplied in non-sterile packaging but sterilized by the user. For before use cleaning and sterilization procedure see chapter "INSTRUCTIONS FOR CLEANING, DRYING AND STERILISATION". The devices are disposable.

G) ANALOGS-Description: The analog is meant to be placed in the impression laboratory model of the mouth of the patient as a reference for the position of the implant in order to create the prosthesis. **Material and Supply:** Manufactured in grade 5 titanium, supplied in non-sterile packaging. The devices are disposable.

H) TRANSFER CAP-Description: The transfer cap is meant to be placed on the impression transfer in order to be incorporated into the impression as a reference for the position of the implant. **Material and Supply:** Manufactured in Peek, supplied in non-sterile packaging but sterilized by the user. For before use cleaning and sterilization procedure see chapter "INSTRUCTIONS FOR CLEANING, DRYING AND STERILISATION". The devices are disposable.

I) CASTABLE ABUTMENT-Description: castable abutments represent a casting element for the creation of dental prostheses in the laboratory. **Materials and supply:** Manufactured from plexiglass, supplied in non-sterile packaging. The devices are disposable.

J) O-RING-Description: The O-Rings allow retention between the o-ball connection of implants or abutments to the metal housing. **Materials and supply:** Nitrile products, supplied in non-sterile packaging. The devices are disposable.

K) SCAN CAP- Description: The scan cap is used on an abutment as a reference point for the implant during intraoral scanning. **Materials and supply:** Made of PEEK. Supplied in non-sterile packaging but sterilized by the user. The device is disposable. For the cleaning and sterilization procedure before use, see the chapter "CLEANING, DRYING AND STERILIZATION INSTRUCTIONS".

LIST OF DRIVERS AND INSTRUMENTS

A) LATCH PROSTHETIC DRIVER-Description: Latch Prosthetic drivers are used to screw or unscrew the prosthetic screw of the prosthesis. **Material and Supply:** Manufactured in stainless steel. Supplied in non-sterile packaging but sterilizable by the user. For before use cleaning and sterilization procedure see chapter "INSTRUCTIONS FOR CLEANING, DRYING AND STERILISATION". The devices are reusable up to a maximum of 30 times. **Torque tightening:** For tightening indications see what is reported in the prosthetic screw section.

B) IMPLANT DRIVER FOR DENTAL IMPLANTS-Description: Implant drivers/mounters are used to place the implant into the surgical site. **Material and Supply:** Manufactured in stainless steel. Supplied in non-sterile packaging but sterilizable by the user. For before use cleaning and sterilization procedure see chapter "INSTRUCTIONS FOR CLEANING, DRYING AND STERILISATION". The devices are reusable up to a maximum of 30 times. **Torque tightening:** 35-40 Ncm for the insertion of the implants both with ratchet and contra angle drivers.

C) IMPLANT MOUNTER- Description: Implant mounters are used to place the implant into the surgical site during guided surgery operations. **Material and Supply:** Manufactured in stainless steel. Supplied in non-sterile packaging but sterilizable by the user. For before use cleaning and sterilization procedure see chapter "INSTRUCTIONS FOR CLEANING, DRYING AND STERILISATION". The devices are reusable up to a maximum of 30 times. **Torque tightening:** 35-40 Ncm for the insertion of the implants both with ratchet and contra angle drivers.

D) LATCH EXTRACTOR- Description: The extractor is used to detach the prosthesis from the EL, CP or BL implants or guided surgery implant mounters from the implant. **Material and Supply:** Manufactured in stainless steel. Supplied in non-sterile packaging but sterilizable by the user. For before use cleaning and sterilization procedure see chapter "INSTRUCTIONS FOR CLEANING, DRYING AND STERILISATION". The devices are reusable up to a maximum of 30 times. **Torque tightening:** Do not exceed 15Ncm for the extraction of the abutment and implant mounters.

Class I Medical Devices

E) MANUAL DRIVERS- Description: Finger Drivers are used to screw or unscrew the prosthetic screw of the prosthesis. **Material and Supply:** Manufactured in stainless steel and Peek. Supplied in non-sterile packaging but sterilizable by the user. For before use and reprocessing procedure see chapter "INSTRUCTIONS FOR CLEANING, DRYING AND STERILISATION". The devices are reusable up to a maximum of 30 times.

F) MANUAL PROSTHETIC DRIVERS- Description: Manual prosthetic drivers are used to screw or unscrew the abutment in position. **Material and Supply:** Manufactured in stainless steel. Supplied in non-sterile packaging but sterilizable by the user. For before use and reprocessing procedure see chapter "INSTRUCTIONS FOR CLEANING, DRYING AND STERILISATION" The devices are reusable up to a maximum of 30 times.

G) DRIVER ADAPTERS-Description: Driver adapters are used to allow a latch driver to be used with a ratchet and a ratchet driver to become a manual driver. **Material and Supply:** Manufactured in stainless steel. Supplied in non-sterile packaging but sterilizable by the user. For before use and reprocessing procedure see chapter "INSTRUCTIONS FOR CLEANING, DRYING AND STERILISATION" The devices are reusable up to a maximum of 30 times.

H) TORQUE RATCHETS-Description: Torque ratchets are used in combination with drivers in order to place in position and tighten dental implants and components. **Material and Supply:** Manufactured in stainless steel and Peek. Supplied in non-sterile packaging but sterilizable by the user. For before use and reprocessing procedure see chapter "INSTRUCTIONS FOR CLEANING, DRYING AND STERILISATION". The devices are reusable up to a maximum of 30 times. **NOTE:** For tightening instructions follow what is indicated in the Implant IFU and in the prosthetic screw section. In order to set the desired torque adjust the nut by rotating it and check the value on the scale. It is advised not to screw beyond the symbol ∞ to avoid stressing internal components excessively, affecting the correct use of the ratchet. For a long service life of the ratchet, after use and sanitation, put it aside with the spring unloaded, namely by adjusting it to the

minimum torque value.

Disassembly instructions (before cleaning process): In order to properly clean the instruments they must be disassembled. The operation can be easily carried out without the use of an additional instrument. No torque ratchet/Peek Ratchet: After use, remove the ratchet wheel and unscrew the driver head so that it can be detached from the handle.

Torque ratchet: pull back the pins in the direction of the body of the ratchet (down) on both sides with your thumb. Then release the pin. Warning: the wheel can be inserted from one side only. Unscrew the torque adjustment nut in order to remove it and disassemble the spring. After cleaning process is over, reassemble the parts and proceed with sterilization.

I) EXTRACTORS - Description: The extractors have the purpose of detaching the prosthetic part from the implant. **Material and supply:** Stainless steel products, supplied in non-sterile packaging but sterilized by the user. The devices are reusable. For the cleaning and sterilization procedure before use, see the chapter "CLEANING, DRYING AND STERILIZATION INSTRUCTIONS".

J) SLEEVES- Description: The sleeve has the purpose of guiding the orientation of the instruments during the surgical operation for greater precision and to fix the surgical guide. **Material and supply:** Manufactured from PEEK and Grade 5 Titanium, supplied in non-sterile packaging. The device is disposable.

K) BOX - Description: The box allows you to collect and organize the medical devices inside. **Material and supply:** Manufactured from polymeric material and aluminium, supplied in non-sterile packaging but sterilized by the user. The devices are reusable. For the cleaning and sterilization procedure before use, see the chapter "CLEANING, DRYING AND STERILIZATION INSTRUCTIONS".

CONTRAINDICATIONS: The contraindication are strictly related to the contraindications of the dental implant system.

General: Uncooperative patient. Intraoral contraindications such as: Anatomical impediments, insufficient bone tissue or poor bone quality. **Related:** Chemotherapy, mild renal disorders, mild liver disorders, endocrine disorders, psychological disorders or psychosis, lack of understanding and motivation, AIDS, seropositive cases, prolonged use of corticosteroids, calcium-phosphorus metabolism disorders, erythropoietic disorders.

Additional contraindications are known allergies to Titanium grade 5, Peek, Stainless Steel, Plexiglass, Chrome Cobalt.

PATIENT SELECTION: Patient selection phase is done during the initial examination of the patient for dental implant suitability.

The medical devices covered by this leaflet are to be used on patients who have been diagnosed with total and/or partial mandibular and/or maxillary edentulism. The device has not been evaluated in paediatric/adolescent patients and is not recommended for use on children. Routine treatment is not recommended until the end of the juvenile growth of the maxillary/mandibular bone has been properly documented. Ethnicity has no impact on the use of these medical devices.

All the raw materials are widely used in the medical field and biocompatible, allergies are rare but possible. It is advisable to check in advance with patients not to have allergies of this type (with particular attention to material Cobalt Chrome, in which case allergies could be more frequent). The manufacturer provides the technical data sheets of the materials upon request for a deeper investigation.

Particular attention should be paid to patients with local or systemic factors that could interfere with or slow down the healing processes of bone or soft tissue and osseointegration (e.g. smoking, poor oral hygiene, uncontrolled diabetes, orofacial radiotherapy, steroid therapy, or infections in the surrounding bone). Special attention should be paid to patients undergoing bisphosphonate therapy. Patients with a hyperfunction (bruxism) may be at greater risk of implant failure or fracture. In general, for patient selection and prosthetic design the dentist should consider all the biological and physical conditions of the patient. The elements for this analysis are acquired during the doctor's studies and career path.

PRE-OPERATIVE PLANNING AND TREATMENT: The suitability of a patient, affected by total or partial edentulism, for implant surgery is derived from a profound analysis of the patient and his clinical history. This analysis is carried out by the surgeon and the dental technician on the basis of the skills acquired during their studies and experience. Once the suitability of the patient for dental implantology has been defined, it is always the task of the surgeon and the dental technician to carry out the surgical and prosthetic planning based on the chosen approach. For further details on the planning of the surgical procedure, consult the information leaflet of the dental implants. The prosthesis that must be made on a single implant or on multiple implants, must be designed using the medical devices supplied by the manufacturer C-Tech Implant in order to fully restore the masticatory and aesthetic functionality of the denture of the patient. Everything must be done in full respect of his/her health and safety.

POST-OPERATIVE RECOVERY: The patient must be informed about the post-operative recovery following surgery as indicated in the dental implant leaflet. Furthermore, the patient must receive all the instructions concerning the correct home oral hygiene practices, the maintenance of the prosthesis (for example removal, in case of removable prostheses, cleaning, disinfection and reassembly) to promote correct patient hygiene and therefore lengthen the life of the device. Regular periodic checkups on hygiene and the state of the prosthesis must be planned by the medical staff.

SIDE EFFECTS: After dental implant surgery, the most common side effects are pain, swelling, speech problems, and gingivitis. Other undesirable effects may result from pre-existing conditions. The dentist should inform the patient of possible common and specific undesirable effects in the case of pre-existing pathologies. Knowledge and identification of specific undesirable effects are a part of a dentist's study programme. The risks of implant surgery include but are not limited to: perforation of the labial or lingual plate, bone fractures, implant fractures, superstructure fractures, aesthetic problems, inadvertent perforation of bone and/or soft tissue (inferior alveolar canal, maxillary sinus, lingual artery, greater palatine foramen), nerve injuries, impairment of the natural dentition, failure or inadequate osseointegration, loss of bone or soft tissue volume due to concomitant events (local infection, peri-implantitis). Pathophysiological problems already diagnosed in the patient may increase the risks described above.

The undesirable effects related to the surgery (surgical phase and healing of the bone and soft tissues) can be consulted in the instruction leaflet of the dental implants. The undesirable effects connected to prosthetic rehabilitation (for example, breakage of the prosthetic screw, breakage of the prosthesis or, in more extreme cases, breakage of the implant) can be linked to inadequate or inaccurate planning, incorrect or improper use of medical device or incorrect positioning of the implants. These undesirable effects can be avoided by following the good practices and knowledge

learned in the course of study, as reported in the leaflets accompanying the manufacturer's medical devices, by consulting the documentation provided and by participating in the courses organized by the manufacturer itself.

DISPOSAL: If it is necessary to dispose of medical devices, particularly when contaminated with blood or fluids, the appropriate hospital waste containers must be used. Use hospital grade gloves to avoid direct contact with the device in question.

INSTRUCTIONS FOR CLEANING, DRYING AND STERILISATION

It is recommended to strictly follow the instructions below for a correct and effective cleaning, disinfection and sterilization procedure. **PLEASE NOTE:** The medical devices supplied by C-Tech Srl as "non-sterile", must be cleaned and sterilised before their first use. The cleaning and sterilisation must be performed according to the methods indicated in the paragraphs below. Reusable devices, before the first use and for all the subsequent uses, must always be subjected to all the cleaning, disinfection, drying and sterilisation steps described below. It is recommended that cleaning operations be carried out as soon as the surgical procedure is completed to minimize the risk of contaminants drying on the devices, making their removal more difficult.

WARNINGS • *Strictly follow the steps described below and carry them out sequentially, without skipping any step.* • *Always use, at all stages, personal protective equipment (such as gloves, mask, goggles, etc.).* • *The devices must not be placed in contact with other objects during the washing phases.* • *Devices that are contained in boxes must be taken out of them and washed separately. The latter must be reinserted in the box, after any reassembly, only after completing all the cleaning phases.* • *Do not use wire brushes or steel wool for cleaning. Residues adhering to instruments, drying up, can cause corrosion.* • *Never use chlorine-based disinfectants.* • *Chemical sterilization is not recommended as this procedure can damage the plastic surface of the surgical boxes and can oxidize the devices.* • *Do not use dry sterilization as the high temperature (about 180 °C) could melt the plastic boxes.* • *Do not sterilize in glass ball sterilizers (temperature above 300 °C).* • *Subject the autoclave, ultrasonic bath and the automated disinfectant to regular maintenance and checks as required by the user manual provided by the manufacturer of these machines.*

CLEANING STEPS- MANUAL OPTION

1. PRE-TREATMENT: Disassemble each part of the multiple instruments (when applicable). Remove the coarsest impurities from the instruments by rinsing with running water for at least 2 minutes. Then remove small residues of blood, saliva, tissue or bone from the instruments using a soft brush with running water, for at least 2 minutes. Do not use brushes with metal or steel wool bristles. The residues adhering to the instruments may cause corrosion as they dry. Such operation must be performed no more than 2 hours after use.

2. CLEANSING: Immediately after rinsing with running water (pre-treatment), immerse disassembled and dirty instruments in a basin with detergent solution for at least one hour, in order to prevent blood residues from coagulating. The instruments must be completely immersed in the solution and the basin must be stirred several times to make the process easier. It is recommended to use a solution based on peracetic acid, such as "SEKUSEPT™ PULVER CLASSIC" with a concentration of 2%. In particular, thoroughly clean the inner side of the cavities of the instruments.

3. DISINFECTION: For a thorough cleaning, place the instruments in the ultrasound machine and leave them immersed for at least 30 minutes and up to an hour. To prevent the corrosion of surgical instruments, such liquids must be suitable for the disinfection of metal and plastic instruments. It is recommended to use a disinfectant based on peracetic acid, such as "SEKUSEPT™ PULVER CLASSIC" with a concentration of 4%. Never use chlorine-based disinfectant. Wash container and instruments separately. Complete the preparation step of steam sterilisation with a thorough and careful rinsing, with sterile water, and then dry as described in the paragraph below.

4. DRYING STEP: Dry the clean instruments in a special hot-air oven at 60 °C for 2 hours. It is recommended to carefully follow the instructions for use provided by the manufacturer of the oven. Make sure that instruments are completely dry, since any residue of steam might start corrosive processes which damage the instrument over time.

AUTOMATED CLEANING & DISINFECTION OPTION

To carry out automated cleaning, disinfection and drying, follow the instructions provided by the manufacturer of the washer disinfectant in use. C-Tech Implant validated the cycle with a washer disinfectant model AT-OS AWD655-10 (DIN 10) that is programmable and equipped with data integrity validated SW. The process is described below.

•Place the devices and any surgical box in the appropriate basket suitable for mechanical washing, close it and load it into the washer disinfectant. •It is recommended to avoid contact of the devices with each other during power washing. •Start the appropriate washing, disinfection and drying programme for the devices. •It requires a pre-rinse with demineralised water for 3 minutes. •Cleaning should be carried out for 5 minutes at a temperature of 55°C, followed by rinsing for 2 minutes at a temperature of 40°C. •We recommend using the concentrated liquid alkaline detergent Neodisher MediClean Forte (0.5%) and the alkaline residue neutraliser Neodisher Z (0.1%). •Final rinse with demineralised water for 2 minutes. •The cycle continues with automated disinfection at 90°C for 5 min. •Take the devices from the washer disinfectant at the end of the drying program and proceed to the sterilisation phase.

STERILISATION STEPS 1) PACKAGING BEFORE STERILISATION -The disinfected, clean and dry instruments must be tidily arranged in the surgical box, which must be sealed in special sterilisable envelopes validated according to ISO 11607-1 and according to what is indicated in the instruction manual supplied by the manufacturer. Alternatively, each instrument can be placed inside paper envelopes (60 g/m²) and plastic film (PET/PP). The manufacturer must indicate that the paper envelopes (60 g/m²) and plastic film (PET/PP) can be steam-sterilised. The paper envelopes and plastic film must be sealed onto all four sides. C-Tech recommends using the same sealing machine model and the same sealing parameters used during the validation phase: LUX SEAL M9025/4 sealing machine; sealing time of 5 seconds; sealing temperature: 150 °C. **2) STERILISATION** For a long service life of the surgical box, sterilisation must be performed in autoclave with steam at a temperature of 121 °C for at least 15 minutes. After sterilization let dry for 20 minutes. It is recommended to validate your autoclave with full and worst-case load in accordance with standard ISO 17665-1. Store the box or the devices individually packaged in a dry environment and inside the relevant envelope for a maximum of 6 months from the sterilisation date. C-Tech Implant recommends using an autoclave model equivalent to the one used during the validation phase: "SERENA 18" model (Manufacturer: Reverberi Srl).

INFORMATION FOR CUSTOMERS: No person is authorised to provide information deviating from the information provided in this instruction sheet. For further information and technical documents, please consult our website at www.c-tech-implant.com.

The Summary of Safety and Clinical Performance (SSCP) is available on EUDAMED, at the following website: <https://ec.europa.eu/tools/eudamed>.

It is connected to the following Basic-UDI:

DIRECT LINK TO INFORMATION MATERIALS: www.c-tech-implant.com

- 805296102PTEMPORARYABTY4 (Plastic Temporary abutment)
- 805296102TTEMPORARYABT5B (Titanium Temporary abutment)
- 805296102FINALABUTMENTQK (Titanium Abutments(straight/angled); Cerec Bases; Dual Channel abutments; Titanium Bases; MUA/FLAT abutments; OMNI abutments; One Time abutments)
- 805296102SPHERICALABT5S (O-Ball Abutments)
- 805296102HEALINGDEVICESND (Healing devices)
- 805296102DENTALTRANSFERCK (Transfer and Scan bodies)
- 805296102DENTALDRIVERSXK (Latch Prosthetic Driver, Implant Driver For Dental Implants, Implant Moulder, Latch Extractor)

If any serious or suspected incidents occur as a result of the use of our medical devices, they should be reported immediately to the Company and to the national competent authority.

“These instructions for use have been translated and are available at the following link: www.c-tech-implant.com. Alternatively, scan the following QRcode”.

In the event of a sudden malfunction of the company website due to force majeure, scan the QRcode to obtain the latest available revision. Alternatively, request the revision of the desired IFU at the following email address quality@c-tech-implant.com. The response will be sent as soon as possible. The paper version of the IFU in your language and which is among those available, can be requested from C-Tech Implant at the email address quality@c-tech-implant.com at any time and will be sent immediately, compatibly with the working hours, and/or within a period not exceeding 7 working days. Alternatively, it can be requested at the time of the order and shipped at the same time, but not inside the packaging. This service is at no additional cost.

SYMBOLS ON THE LABELS OF MEDICAL DEVICES

Table 1: Definition of the symbols on the labels.

Symbol	Description
	CE Marking of a Medical Device If applicable: the identification number of the notified body must follow this symbol.
	CE marking of a Medical Device with intervention of the Notified Body TÜV Rheinland
	Code of the device
	Lot Number
	Manufacturer
	Manufacturing date
	Medical device sterilized using irradiation. Single sterile barrier system with protective packaging outside
	Non Sterile device
	Disposable device
	Do not re-sterilise
	Expiry date
	Warning
	Read the package leaflet carefully or consult the electronic IFU here: www.c-tech-implant.com
	Do not use if packaging is damaged or opened


	Keep away from direct sunlight
	Keep in a dry place
	Medical device
	Contains Hazardous substances
	Unique Device Identifier

Table 2

Symbol	Description
	Patient Identification
	Information website for patient
	Nome of the doctor or the healthcare centre
	Date of the surgery
	Medical Device

UDI Code and Data Matrix

(01): UDI-DI, identifies the manufacturer and the specific version or model of device. (10): Batch number. (17): Expiry date. (11): Date of production.

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